

National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK)

Polypeptide Hormones, Subunits, and Antisera Available to Qualified Investigators for Research

The materials listed below are procured or donated through programs supported by the National Institute of Diabetes and Digestive and Kidney Diseases, the Center for Population Research of the National Institute of Child Health and Human Development, the Agricultural Research

Service of the U.S. Department of Agriculture and Dr. A. F. Parlow of the Harbor-UCLA Medical Center. These materials are distributed to qualified investigators for research purposes, not for therapeutic, diagnostic or commercial uses.

Application—Application forms must be submitted prior to an award. Letters or requests not accompanied by an application form are not acceptable. Forms are available from and should be submitted to:

Terry Taylor, MD, Director, National Hormone and Pituitary Program
685 Lofstrand Lane, Rockville, MD 20850
Voice PHONE: (301) 309-3667, FAX PHONE: (301) 340-9245

Awards—Only one ampoule or one RIA kit of each material will be awarded at a time. Requests for larger amounts of materials may be

granted, subject to review by the Hormone Distribution Program Advisory Committee.

New Materials Available:

- Human pituitary FSH, LH, and TSH: Free of hormonal contaminants
- New, monospecific polyclonal antisera to rat PRL, GH, FSH β , LH β and TSH β
- New, sensitive antiserum for rat LH RIA
- Rat GH and PRL free of hormonal contaminants
- Ovine, bovine, porcine FSH devoid of LH
- Rat LH, FSH and TSH subunits/antisera
- Glycosylated ovine and porcine PRL
- 20K variant of hGH
- Bovine FSH RIA
- Antiserum to cAMP

For RIA: Pituitary and placental hormones and subunits, antisera and reference preparations

	Human			Rat			Ovine			Bovine		Porcine			Monkey		
	H	AS	RP	H	AS	RP	H	AS	RP	H	AS	H	AS	RP	H	AS	RP
ACTH	1a	1	—	—	—	1	—	—	—	—	—	—	—	—	—	—	—
CG	Intact	1b	1	—	—	—	—	—	—	—	—	—	—	—	—	—	—
CG α	Subunit	1a	1	—	—	—	—	—	—	—	—	—	—	—	—	—	—
CG β	Subunit	1b	1	—	—	—	—	—	—	—	—	—	—	—	—	—	—
FSH	Intact	1b	1	1	1b	1	1	1a	1	1	1	1o	1c	1p	—	1b	1
FSH α	Subunit	1a	1	—	1b	1	1	*	—	—	—	—	*	*	—	—	—
FSH β	Subunit	1b	1	—	1b	1j	1	*	—	—	—	—	*	*	—	—	—
GH		1h	1	1	1e	1	1	1h	1	—	*	*	*	*	—	*	—
GH, 20K	Variant	1c,g	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—
LH	Intact	1e	1	1	1b	1	1	1f	1	—	1g	—	1c	*	—	1n	1
LH α	Subunit	1a	1k	—	1b	1	1	1c	1	—	1	—	*	—	—	—	—
LH β	Subunit	1c	1	—	1b	1	1	1c	1m	—	1	—	*	*	—	—	—
LPH β		1d	1	1	—	—	—	—	—	—	—	—	—	—	—	—	—
PL		1i	1	—	—	—	—	—	—	—	—	—	—	—	—	—	—
PRL		1b	1	1	1f	1	1	1g	1	—	1g	*	1b	*	—	—	—
TSH	Intact	1e	1	1	1b	1	1	1g	1	—	1b	*	*	—	—	—	—
TSH α	Subunit	1c	1	—	1b	1	1	*	—	—	*	—	*	—	—	—	—
TSH β	Subunit	1b	1l	—	1b	1	1	*	—	—	1	*	*	—	—	—	—
(PTH)		—	—	—	—	—	—	—	—	—	2 μ g	—	—	—	—	—	—

a 50 μ g; b 100 μ g; c 200 μ g; d 300 μ g; e 400 μ g; f 500 μ g; g 1mg; h 2mg; i 5mg; j Antisera for 2,500 RIAs; k Antisera for 3,750 RIAs; l Antisera for 5,000 RIAs; m Antisera for 9,750 RIAs; n These materials will be shipped frozen by air express in the U.S.A.; o Request antiserum to ovine FSH, for bovine FSH RIA; p Request antiserum to ovine FSH, for porcine FSH RIA. * These may be obtained from Dr. A. F. Parlow, Harbor-UCLA Medical Center, 1000 West Carson Street, Torrance, CA 90509. H, Hormone for Iodination; AS, Diluted Antiserum for 10,000 RIAs or as indicated; RP, Reference Preparation.

Note—RIAs for ovine pituitary hormones cross-react strongly with, and may be used for the measurement of, the corresponding bovine pituitary hormones, with the exception of FSH. RIAs for rat pituitary hormones cross-react strongly with, and may be used for the measurement of, the corresponding mouse pituitary hormones, with the exception of PRL. For mouse PRL, use the available homologous

RIA. For mouse GH RIA, use the available mouse GH for iodination with rat GH antiserum. RIAs for human pituitary GH, PRL, and TSH cross react strongly with the corresponding monkey pituitary hormones. Antisera are distributed in liquid or lyophilized form. Antisera for RIA contain merthiolate, if distributed in liquid form.

For biologic studies: Polypeptide hormones and other materials

	Human	Rat	Ovine	Bovine	Porcine	Other
CG	1 mg	—	—	—	—	—
CGa	1 mg	—	—	—	—	—
CG hybrid (ABIER-CR-1xy)	1.5 mg	—	—	—	—	—
FSH	2 mg	—	0.50 mg	—	10 mg	—
GH	10 mg	5 mg	50 mg	10 mg	10 mg	—
GH, 20K Variant	1 mg	—	—	—	—	—
LH	5 mg	—	10 mg	5 mg	10 mg	—
PRL	2 mg	5 mg	50 mg	10 mg	5 mg	—
PRL, Glycosylated	—	—	1 mg	—	1 mg	—
TSH	2 mg	—	1 mg	1 mg	—	—
Folliculostatin Reference Preparation <i>n</i>	—	—	—	—	0.2 U	—
Insulin (guinea pig)	—	—	—	—	—	80 µg
Insulin (desoctapeptide)	—	—	—	—	100 µg	—
LH-FSH-RH (GnRH)	—	—	—	—	—	100 µg
LH-FSH-RH (GnRH)	—	—	—	—	—	10 mg
Porcine follicular fluid (bottle) <i>n</i>	—	—	—	—	500 ml	—
PTH (BIO)	—	—	—	1 mg	—	—
Relaxin	—	—	—	—	1 mg	—
Follistatin (recombinant)	0.1 mg	—	—	—	—	—
Follistatin (recombinant)	1 mg	—	—	—	—	—

PITUITARY GLANDS, FROZEN, HUMAN

U.S. investigators only. Not for clinical investigation.

ANTISERA FOR IMMUNOCYTOCHEMISTRY (undiluted or in low dilution)

hGH natural, hGH synthetic, hPRL, hFSH β , hLH β , hTSH β , and hACTH	1 ampoule
rGH, rPRL, rFSH β , rLH β , rTSH β and mouse GH	1 ampoule
rEstrogen Receptor (ER) Antibody and Synthetic ER Peptide (amino acids 270–284)	1 ampoule each
oLH β and oPRL	1 ampoule

OTHER ANTISERA

Aldosterone antiserum	1 ampoule
Somatomedin-C antiserum (polyclonal)	1 ampoule
Somatomedin-C antiserum (monoclonal)	1 ampoule

HYPOTHALAMI, LYOPHILIZED, RAT

MATERIALS AVAILABLE THROUGH THE COURTESY OF DR. A. F. PARLOW

The antisera listed below may be obtained by writing directly to Dr. A. F. Parlow, Director, Pituitary Hormones and Antisera Center, Harbor–UCLA Medical Center, 1000 West Carson Street, Torrance, CA 90509.

Phone: (310) 222-3537; FAX Phone: (310) 222-3432; BITNET: "PARLOW@HARBORUC.BITNET"

RIA Kits	Canine	Chicken	Equine	Hamster	Mouse	Rabbit
GH	1	1	1	—	1	1
PRL	1	1	1	1	1	1
FSH	1	—	—	—	RP	1
LH	1	—	RP	—	1	1
TSH	—	—	1	—	RP	1

1, 1 kit; RP, Reference Preparation.

FOR MATERIALS OTHER THAN THOSE LISTED ABOVE, INQUIRE OF DR. TERRY TAYLOR OR DR. A. F. PARLOW.

Taylor & Francis

Journals of quality since 1798

ARCHIVES OF ANDROLOGY

Edited by:

*E. S. E. Hafez
Reproductive Health Center, USA*

ARCHIVES OF ANDROLOGY publishes research papers, review articles, and short communications in the areas of reproduction, regulation and infertility in the human and animal male. Aspects covered include: endocrinological, morphological, microbiological, neurological, anatomical, physiological, immunological, ultrastructural, pathological, psychiatric, fertility regulation and contraception, micromanipulation/cryopreservation of gametes, HIV / STD (sexually transmitted disease), in vitro fertilization / embryo transfer (IVF / ET).

Readership

Physiologists, andrologists, pathologists, endocrinologists, immunologists, psychiatrists, neurologists, obstetricians, gynaecologists and fertility specialists.

Subscription Information

Vol 30/31 (1993), Bimonthly
ISSN: 0148-5016
Institutional: US\$405 / £233
Personal: US\$180 / £128

**Send today for your FREE sample
copy!**

UK:

Rankine Road, Basingstoke,
Hampshire RG24 0PR

Tel: (0256) 840366 Fax: (0256) 479438

USA:

1900 Frost Road, Suite 101,
Bristol, PA 19007-1598


Tel: (215) 785 5800 Fax: (215) 785 5515



Taylor & Francis
Publishers since 1798

MAKLER® COUNTING CHAMBER

For Rapid Sperm Analysis From Undiluted Specimen



Sefi-Medical Instruments Ltd.

P.O.B. 7295, Haifa 31070, Israel, Tel. (972)-4-251651
Telex: 46400 IL. Ext. 8796, Fax: (972)-4-258903

ADVANTAGES

- **SPERM ANALYSIS PERFORMED FROM UNDILUTED SPECIMEN**
- **EASY AND FAST RESULTS:** The number of spermatozoa counted in any strip of 10 squares of the grid indicates their concentration as millions/ml. **No factors necessary for calculations.**
- **OPTIMAL DEPTH:** 10 microns. **No blurring.** Free lateral head movement. Uniform mono-layer in one focal plane.
- **BUILT-IN GRID:** Grid on cover glass. No insertion of grid into microscope eye piece.
- **ECONOMICAL: Reusable.** Easily wiped with non-bleach disinfectant solution.
- **SPERM COUNT, MOTILITY EVALUATION AND MORPHOLOGY DETERMINATION ALL IN ONE CHAMBER**
- **CALIBRATION UNNECESSARY:** **Repeated use** with no change in accuracy due to quartz-coated contact points.
- **SUPERIOR ACCURACY:** State-of-the-art precision engineering. Laser beam checked.

THE MAKLER® INSEMINATION DEVICE

FOR INTRA-UTERINE INSEMINATION



This easy-to-use device with a disposable atraumatic cannula, has been specially designed to prevent leakage of the injected contents.

Applications:

- ★ For infertility due to subfertile oligo-asthenospermic specimens.
- ★ For cervical mucus hostility, low semen volume or anti-sperm antibodies.
- ★ In conjunction with induction of ovulation, to increase probability of fertilization.
- ★ In case of unexplained infertility.
- ★ For increased efficacy in donor insemination, or the use of one specimen for more than one patient.




Sefi-Medical Instruments

P.O.B. 7295, HAIFA 31070, ISRAEL. TEL. (972)-4-251651
TELEX: 46400 IL. EXT. 8796. FAX. (972)-4-258903

New HEIGHTS in Flexible cell analysis

Building upon the rock solid reliability of the **HTM-IVOS**, the new Version 10 system climbs to innovative heights in semen analysis. With multiple magnifications, extended analysis options, and new Windows™ based software, the Version 10 HTM-IVOS is the pinnacle of versatility.



HTM-IDENT

HTM-IDENT, the first option of the Spectra series of integrated, fluorescent technologies, yields a higher degree of accuracy in sperm counts and concentrations. With the use of a specialized DNA-specific dye, the sperm are readily identified under fluorescent illumination and easily discriminated from debris.

HTM-DIMENSIONS

Expanding into the field of morphology analysis, **HTM-DIMENSIONS** automated analysis system objectively classifies sperm by strict analysis of the head, acrosome, and midpiece. High power, 100x magnification produces crisp images of stained sperm for analysis.

Call today for more information about the new Version 10 HTM-IVOS, HTM-IDENT, and HTM-DIMENSIONS, the three new "Innovations To Rely On" from Hamilton-Thorn Research.