


Verreault
MARCH/APRIL 1993
Volume 14 Number 2

Journal of
ANDROLOGY



Published by THE AMERICAN SOCIETY OF ANDROLOGY

ACCUDENZ®

DENSITY GRADIENT MEDIUM

For The Selection of Motile Human Spermatozoa

- Highest % of motility up to 21 hours after selection
- Better selection of motile spermatozoa from abnormal semen
- Out-performs the swim-up method in the selection of motile spermatozoa in male factor patients (see chart)
- Non-toxic

References: "Sperm Selection Using Density Gradients" Fertility & Sterility, February 1988 issue, page 335.


ACCUDENZ®

DENSITY GRADIENT MEDIUM

For Study of Buoyancy and Sedimentation of Human X- and Y-bearing sperm

References: Archives of Andrology 19:211-215 (1987)

MOTILITY TEST RESULTS:



% of Motile Spermatozoa
21 Hours After Selection From
Male Factor Patients Using Three
Different Methods.

Accudenz®...
THE PERFORMANCE MEDIUM.

ACCUDENZ®

The Autoclavable, Universal Centrifugation Medium!

Ordering #	Description	Packing	Price
AN 7050	Accudenz A.G. (Analytical Grade)	1 x 50 gm. 10 x 50 gm.	Send For New Low Prices
AN 7050/BLK	Accudenz A.G.	1 x 500 gm.	

**NEW
LOW
PRICES!**

Send for Cell Separation Media brochure and applications chart.
Accudenz® technical booklets and FREE samples available on request.

For research purposes only. Not for human or diagnostic use.

© Copyright 1993 Accurate Chemical & Scientific Corp.®

ACCURATE Chemical & Scientific Corporation

300 SHAMES DRIVE, WESTBURY, NY 11590 USA • TELEX: 4972582

FAX: 1-516-997-4948

WESTBURY, NY: 516-333-2221

SAN DIEGO, CA: 619-235-9400


TOLL FREE: 1-800-645-6264

OR

TOLL FREE: 1-800-255-9378

MAKLER® COUNTING CHAMBER

For Rapid Sperm Analysis From Undiluted Specimen



Sefi-Medical Instruments Ltd.


P.O.B. 7295, Haifa 31070, Israel, Tel. (972)-4-251651
Telex: 46400 IL. Ext. 8796, Fax: (972)-4-258903

ADVANTAGES

- **SPERM ANALYSIS PERFORMED FROM UNDILUTED SPECIMEN**
- **EASY AND FAST RESULTS:** The number of spermatozoa counted in any strip of 10 squares of the grid indicates their concentration as millions/ml. **No factors necessary for calculations.**
- **OPTIMAL DEPTH:** 10 microns. **No blurring.** Free lateral head movement. Uniform mono-layer in one focal plane.
- **BUILT-IN GRID:** Grid on cover glass. No insertion of grid into microscope eye piece.
- **ECONOMICAL: Reusable.** Easily wiped with non-bleach disinfectant solution.
- **SPERM COUNT, MOTILITY EVALUATION AND MORPHOLOGY DETERMINATION ALL IN ONE CHAMBER**
- **CALIBRATION UNNECESSARY:** **Repeated use** with no change in accuracy due to quartz-coated contact points.
- **SUPERIOR ACCURACY:** State-of-the-art precision engineering. Laser beam checked.

THE MAKLER® INSEMINATION DEVICE

FOR INTRA-UTERINE INSEMINATION



This easy-to-use device with a disposable atraumatic cannula, has been specially designed to prevent leakage of the injected contents.


Applications:

- ★ For infertility due to subfertile oligo-asthenospermic specimens.
- ★ For cervical mucus hostility, low semen volume or anti-sperm antibodies.
- ★ In conjunction with induction of ovulation, to increase probability of fertilization.
- ★ In case of unexplained infertility.
- ★ For increased efficacy in donor insemination, or the use of one specimen for more than one patient.



Sefi-Medical Instruments

P.O.B. 7295, HAIFA 31070, ISRAEL. TEL. (972)-4-251651
TELEX: 46400 IL. EXT. 8796. FAX. (972)-4-258903



Setting the Standard for Excellence

Quality First!

- When you want reliability
- When you want consistency
- When you want convenience

...In Sperm Preparation



IRVINE SCIENTIFIC

Santa Ana, CA

For Further Information
or Customer Services

Contact:

714-261-7800

800-437-5706

FAX: 714-261-6522

Journal of ANDROLOGY

Professional Position Announcements

You won't find a more economical and effective way to reach qualified candidates for professional positions in the field of andrology. Position listings may be submitted as 1/4, 1/2 and full page camera-ready ads. Typesetting is available for an additional charge. For rates and deadlines, please contact:

Journal of Andrology Advertising Office
1041 New Hampshire Street
Lawrence, KS 66044 USA
Tel: (800) 627-0932 • Fax: (913) 843-1274

New HEIGHTS in Flexible cell analysis

Building upon the rock solid reliability of the **HTM-IVOS**, the new Version 10 system climbs to innovative heights in semen analysis. With multiple magnifications, extended analysis options, and new Windows™ based software, the Version 10 **HTM-IVOS** is the pinnacle of versatility.



HTM-IDENT

HTM-IDENT, the first option of the Spectra series of integrated, fluorescent technologies, yields a higher degree of accuracy in sperm counts and concentrations. With the use of a specialized DNA-specific dye, the sperm are readily identified under fluorescent illumination and easily discriminated from debris.

HTM-DIMENSIONS

Expanding into the field of morphology analysis, **HTM-DIMENSIONS** automated analysis system objectively classifies sperm by strict analysis of the head, acrosome, and midpiece. High power, 100x magnification produces crisp images of stained sperm for analysis.

Call today for more information about the new Version 10 HTM-IVOS, HTM-IDENT, and HTM-DIMENSIONS, the three new "Innovations To Rely On" from Hamilton-Thorn Research.

Quick, easy measurement of testes.


An orchidometer that is accepted around the world comes to the USA!

The measurement of testicular size is greatly simplified with this new orchidometer which has been in use worldwide for over 25 years and is now made available in America by ASSI®.

It consists of a graduated series of 12 models, each precisely defined in shape and volume, that allows the pediatrician, endocrinologist or andrologist to make the necessary comparison and evaluation.

This orchidometer is both simple to use and safe, completely eliminating the use of calipers.

Call or write for full information.



**accurate
surgical &
scientific
instruments
corporation**

300 SHAMES DRIVE, WESTBURY, NY 11590
TELEPHONE: (516) 333-2570
FAX: 1-516-997-4948

TO ORDER CALL TOLL FREE:

1-800-645-3569

OR

1-800-255-9378

(WEST COAST OFFICE)

ASSI®

© Copyright 1993 ASSI®

ACCURATE SURGICAL & SCIENTIFIC
INSTRUMENTS CORPORATION
300 Shames Drive, Westbury, NY 11590

Please send me more information on the orchidometer.

Name _____

Institution _____

Address _____

City, State, Zip _____

Phone No. _____

ANDR