



MAY/JUNE 1995  
Volume 16 Number 3

# Journal of ANDROLOGY



Published by THE AMERICAN SOCIETY OF ANDROLOGY

# ACCUDENZ<sup>®</sup>

## DENSITY GRADIENT MEDIUM

### For The Selection of Motile Human Spermatozoa

- Highest % of motility up to 21 hours after selection
- Better selection of motile spermatozoa from abnormal semen
- Out-performs the swim-up method in the selection of motile spermatozoa in male factor patients (see chart)
- Non-toxic


References: "Sperm Selection Using Density Gradients" Fertility & Sterility, February 1988 issue, page 335.

# ACCUDENZ<sup>®</sup>

## DENSITY GRADIENT MEDIUM

### For Study of Buoyancy and Sedimentation of Human X- and Y-bearing sperm

References: Archives of Andrology 19:211-215 (1987)



% of Motile Spermatozoa  
21 Hours After Selection From  
Male Factor Patients Using Three  
Different Methods.

Accudenz<sup>®</sup>...  
**THE PERFORMANCE MEDIUM.**

# ACCUDENZ<sup>®</sup>

## The Autoclavable, Universal Centrifugation Medium!

Ordering #	Description	Packing	Price
AN 7050	Accudenz A.G. (Analytical Grade)	1 x 50 gm. 10 x 50 gm.	Send For New Low Prices
AN 7050/BLK	Accudenz A.G.	1 x 500 gm.	

Send for Cell Separation Media brochure and applications chart.  
Accudenz<sup>®</sup> technical booklets and FREE samples available on request.

For research purposes only. Not for human or diagnostic use.

© Copyright 1995 Accurate Chemical & Scientific Corp.

**ACCURATE** Chemical & Scientific Corporation

300 SHAMES DRIVE, WESTBURY, NY 11590 USA • FAX: 1-516-997-4948


WESTBURY, NY: 516-333-2221 OR SAN DIEGO, CA: 619-296-9945

TOLL FREE: 1-800-645-6264 OR TOLL FREE: 1-800-255-9378

®

# MAKLER® COUNTING CHAMBER

For Rapid Sperm Analysis From Undiluted Specimen



**Sefi-Medical Instruments Ltd.**

P.O.B. 7295, Haifa 31070, Israel, Tel. (972)-4-251651  
Telex: 46400 IL. Ext. 8796, Fax: (972)-4-258903

## ADVANTAGES

- **SPERM ANALYSIS PERFORMED FROM UNDILUTED SPECIMEN**
- **EASY AND FAST RESULTS:** The number of spermatozoa counted in any strip of 10 squares of the grid indicates their concentration as millions/ml. **No factors necessary for calculations.**
- **OPTIMAL DEPTH:** 10 microns. **No blurring.** Free lateral head movement. Uniform mono-layer in one focal plane.
- **BUILT-IN GRID:** Grid on cover glass. No insertion of grid into microscope eye piece.
- **ECONOMICAL: Reusable.** Easily wiped with non-bleach disinfectant solution.
- **SPERM COUNT, MOTILITY EVALUATION AND MORPHOLOGY DETERMINATION ALL IN ONE CHAMBER**
- **CALIBRATION UNNECESSARY:** **Repeated use** with no change in accuracy due to quartz-coated contact points.
- **SUPERIOR ACCURACY:** State-of-the-art precision engineering. Laser beam checked.

# THE MAKLER® INSEMINATION DEVICE

FOR INTRA-UTERINE INSEMINATION



This easy-to-use device with a disposable atraumatic cannula, has been specially designed to prevent leakage of the injected contents.

#### Applications:

- ★ For infertility due to subfertile oligo-asthenospermic specimens.
- ★ For cervical mucus hostility, low semen volume or anti-sperm antibodies.
- ★ In conjunction with induction of ovulation, to increase probability of fertilization.
- ★ In case of unexplained infertility.
- ★ For increased efficacy in donor insemination, or the use of one specimen for more than one patient.



**Sefi-Medical Instruments**

P.O.B. 7295, HAIFA 31070, ISRAEL. TEL. (972)-4-251651  
TELEX: 46400 IL. EXT. 8796. FAX. (972)-4-258903

***Free On-Line Document Search...  
Discount on Document Delivery***



**UnCover**

**Announcing  
*Access to UnCover***

**Document search and delivery for readers  
of *Journal of Andrology***

You can now access free of charge the UnCover on-line database of the complete table of contents information from over 16,000 journals. In addition, you can order delivery of full text of articles at a discount of \$2.00 per article off the regular document delivery price—a discount normally available only to UnCover subscribers who pay an annual password fee.

With *Access to UnCover*, you can search for articles on particular topics, by particular authors, or from particular journals. You can then order on-line the full text of articles for delivery by fax within 24 hours. Document delivery fees will be billed to your credit card.

*Access to UnCover* is available now over the Internet.

Simply Telnet to the *Access to UnCover* Address:

allen.carl.org

At the prompt, enter the Password: DELIVER

***Access to UnCover* makes it easy to find  
the articles you need.**


# Quick, easy measurement of testes.

Now available from ASSI®, the Test-Size® orchidometer that is accepted around the world!

The measurement of testicular size is greatly simplified with this new orchidometer which has been in use worldwide for over 25 years and is now

made available in America by ASSI®.


It consists of a graduated series of 12 models, each precisely defined in shape and volume, that allows the pediatrician, endocrinologist or andrologist to make the necessary comparison and evaluation.



This orchidometer is both simple to use and safe, completely eliminating the use of calipers.

**Sterilization information: Steam autoclave at 250°F for 15 to 20 minutes.**

Test-Size® is a registered trademark of Remcat Trade A.B.  
© Copyright 1993 ASSI®



**accurate surgical & scientific instruments corporation**

300 SHAMES DRIVE, WESTBURY, NY 11590  
TELEPHONE: (516) 333-2570  
FAX: 1-516-997-4948

TO ORDER CALL TOLL FREE:

1-800-645-3569

OR

1-800-255-9378

(WEST COAST OFFICE)

*Exclusive North American distributor of the Test-Size® orchidometer.*

**ASSI®**  
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS™

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS CORPORATION  
300 Shames Drive, Westbury, NY 11590

Please send me more information on the Test-Size® orchidometer.

Name \_\_\_\_\_

Institution \_\_\_\_\_

Address \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Phone No. \_\_\_\_\_

ANDR