



SEPTEMBER/OCTOBER 1995  
Volume 16 Number 5

Journal of  
**ANDROLOGY**



Published by THE AMERICAN SOCIETY OF ANDROLOGY

© Copyright 1995 Accurate Chemical & Scientific Corp.




# YES, WE HAVE MONOCLONALS



Look for our products in Linscott's Directory

## WE HAVE <sup>NEW</sup> MONOCLONALS TODAY!

**ACCURATE Chemical & Scientific Corporation**  
300 SHAMES DRIVE, WESTBURY, NY 11590 USA • FAX: 1-516-997-4948  
WESTBURY, NY: 516-333-2221    SAN DIEGO, CA: 619-296-9945  
TOLL FREE: 1-800-645-6264    OR    TOLL FREE: 1-800-255-9378




ACCURATE CHEMICAL & SCIENTIFIC CORPORATION  
300 SHAMES DRIVE, WESTBURY, N.Y. 11590 USA  
 Please send me the NEW "Incomplete" Accurate catalog.  
NAME \_\_\_\_\_  
INSTITUTION \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
CITY, STATE & ZIP \_\_\_\_\_

®

ANDR

# New HEIGHTS in Flexible cell analysis

Building upon the rock solid reliability of the **HTM-IVOS**, the new Version 10 system climbs to innovative heights in semen analysis. With multiple magnifications, extended analysis options, and new Windows™ based software, the Version 10 **HTM-IVOS** is the pinnacle of versatility.



## **HTM-IDENT**

**HTM-IDENT**, the first option of the Spectra series of integrated, fluorescent technologies, yields a higher degree of accuracy in sperm counts and concentrations. With the use of a specialized DNA-specific dye, the sperm are readily identified under fluorescent illumination and easily discriminated from debris.


## **HTM-DIMENSIONS**

Expanding into the field of morphology analysis, **HTM-DIMENSIONS** automated analysis system objectively classifies sperm by strict analysis of the head, acrosome, and midpiece. High power, 100x magnification produces crisp images of stained sperm for analysis.

**Call today for more information about the new Version 10 HTM-IVOS, HTM-IDENT, and HTM-DIMENSIONS, the three new "Innovations To Rely On" from Hamilton-Thorn Research.**

# MAKLER® COUNTING CHAMBER

For Rapid Sperm Analysis From Undiluted Specimen



**Sefi-Medical Instruments Ltd.**

P.O.B. 7295, Haifa 31070, Israel, Tel. (972)-4-251651  
Telex: 46400 IL. Ext. 8796, Fax: (972)-4-258903

## ADVANTAGES

- **SPERM ANALYSIS PERFORMED FROM UNDILUTED SPECIMEN**
- **EASY AND FAST RESULTS:** The number of spermatozoa counted in any strip of 10 squares of the grid indicates their concentration as millions/ml. **No factors necessary for calculations.**
- **OPTIMAL DEPTH:** 10 microns. **No blurring.** Free lateral head movement. Uniform mono-layer in one focal plane.
- **BUILT-IN GRID:** Grid on cover glass. No insertion of grid into microscope eye piece.
- **ECONOMICAL: Reusable.** Easily wiped with non-bleach disinfectant solution.
- **SPERM COUNT, MOTILITY EVALUATION AND MORPHOLOGY DETERMINATION ALL IN ONE CHAMBER**
- **CALIBRATION UNNECESSARY:** **Repeated use** with no change in accuracy due to quartz-coated contact points.
- **SUPERIOR ACCURACY:** State-of-the-art precision engineering. Laser beam checked.

# THE MAKLER® INSEMINATION DEVICE

FOR INTRA-UTERINE INSEMINATION



This easy-to-use device with a disposable atraumatic cannula, has been specially designed to prevent leakage of the injected contents.

#### **Applications:**

- ★ For infertility due to subfertile oligo-asthenospermic specimens.
- ★ For cervical mucus hostility, low semen volume or anti-sperm antibodies.
- ★ In conjunction with induction of ovulation, to increase probability of fertilization.
- ★ In case of unexplained infertility.
- ★ For increased efficacy in donor insemination, or the use of one specimen for more than one patient.



**Sefi-Medical Instruments**

P.O.B. 7295, HAIFA 31070, ISRAEL. TEL. (972)-4-251651  
TELEX: 46400 IL. EXT. 8796. FAX. (972)-4-258903

# Quick, easy measurement of testes.

*Now available from ASSI®, the Test-Size® orchidometer that is accepted around the world!*

*The measurement of testicular size is greatly simplified with this new orchidometer which has been in use worldwide for over 25 years and is now*


*made available in America by ASSI®.*

*It consists of a graduated series of 12 models, each precisely defined in shape and volume, that allows the pediatrician, endocrinologist or andrologist to make the necessary comparison and evaluation.*

*This orchidometer is both simple to use and safe, completely eliminating the use of calipers.*

*Sterilization information: Steam autoclave at 250°F for 15 to 20 minutes.*

Test-Size® is a registered trademark of Remcat Trade A.B.  
© Copyright 1993 ASSI®



**accurate surgical & scientific instruments corporation**

300 SHAMES DRIVE, WESTBURY, NY 11590  
TELEPHONE: (516) 333-2570  
FAX: 1-516-997-4948

TO ORDER CALL TOLL FREE:

**1-800-645-3569**

OR

**1-800-255-9378**

(WEST COAST OFFICE)

*Exclusive North American distributor of the Test-Size® orchidometer.*

**ASSI®**  
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS™

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS CORPORATION  
300 Shames Drive, Westbury, NY 11590

Please send me more information on the Test-Size® orchidometer.

Name \_\_\_\_\_

Institution \_\_\_\_\_

Address \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Phone No. \_\_\_\_\_

ANDR